



Made in Italy

OVER MASK-I

Reusable Filtering Half Mask



1. Wash hands and face with soap before and after using the mask and before touching it




2. Wash mask before use and before wearing it





3. Take the mask from the external part, put the mask on the face-covering your nose, you mouth and your chin: be sure the elastics are well places on your ears.



4. After putting the mask properly, do not touch the mask with your hands.

Size	Single size
Code	M035-B100 BOC containing 100 pcs
	M035-K100 CARTON containing 1000 pcs. (10 boxes containing 100 pcs)
Single Packed Mask	 The single packaging improves hygiene by avoiding any contact with harmful substances before the use.

MATERIALS

OUTER LAYER	 100% water repellent polypropylene (70g/m ² ± 10%) Tested for harmful substances
INNER LAYER	 100% water repellent polypropylene (70g/m ² ± 10%) Tested for harmful substances
ELASTIC	Spandex and nylon

TECHNICAL SPECIFICATIONS

TESTS	RESULTS
BREATHABILITY WITH WATER VAPOR	20,4 mg/cm ² *h
WATERPROOFNESS	0 mm water dragging



MAXIMUM TEMPERATURE 60 C (THREE WASHING CYCLES ARE GUARANTEED)

FOR A LONGER LIFE OF THE PRODUCT, GENTLY HAND WASH AT 35-40 C



DO NOT IRON

